

PDV at 532 nm

PDV Workshop, May 2018

Patrick Younk
LA-UR 18-24003



Operated by Los Alamos National Security, LLC for the U.S. Department of Energy's NNSA

The Plan of the Talk

- Why is using a different frequency for PDV interesting?
- Prototype System and Early Results
- Challenges
- Conclusions

Note: This work has been supported by funding from Laboratory Directed Research and Development at Los Alamos National Lab.

Why is using a different frequency for PDV interesting?

At least four independent reasons:

1. The resolution of PDV is proportional to λ .
2. The scatter from particles is a strong function of the size parameter D/λ .
3. The diffraction limited spot size is proportional to λ .
4. Higher frequency light can penetrate denser plasmas.


I will review reasons 1 and 2 here.

Confirming reasons 3 and 4 are homework for you.

Possible Applications for Shorter Wavelength PDV

1. Detonator Function Time
2. Grain Selective Measurements
3. Coherent Detection of sub-micron particles (ejecta)
4. Warm Dense Matter

Recent Paper: Chaoqun Jiang, et al., “A 532 nm fiber-optic displacement interferometer for low-velocity impact experiments”, Review of Scientific Instruments 89, 023101 (2018).




This paper was published out about 18 months after I submitted my proposal

The resolution of PDV is proportional to λ .

A bit of a review:

Resolution is the ability to separate two different events (time resolution) or two different surfaces (velocity or position resolution).



Resolution of Motion is proportional to λ .

Typical Resolution of Motion:

For many of my 1550 nm PDV experiments we quote:

$$\Delta t = 12 \text{ ns} \text{ and } \Delta v = 12 \text{ m/s}$$

Checking with Equation #1:

$$\Delta t \Delta v \approx \lambda / (4\pi)$$

$$(12 \text{ ns}) (12 \text{ nm/ns}) \approx 1550 \text{ nm} / (4\pi) \checkmark$$

For PDV at 532 nm we could expect:


$$\Delta t = 3 \text{ ns} \text{ and } \Delta v = 12 \text{ m/s} \text{ or}$$

$$\Delta t = 12 \text{ ns} \text{ and } \Delta v = 3 \text{ m/s, for example.}$$

Particle detection with PDV....

- **First Point: PDV is a coherent detection method and has excellent performance with small signal returns.**
- **Second Point:**
 - The cross Section of a particle in the Rayleigh Regime is proportional to D^6/λ^4
 - The strong dependence on λ , means that you get a big gain when going to shorter wavelengths.
 - The strong dependence on D , suggests a natural way to constrain the distribution of particles sizes.
- **Trivia: Weather Radar is coherent detection of ~1 mm particles (rain drops) with 10 cm light, which can quantify areal mass.**


Prototype System



Benchtop Setup



Probe and Target



Early Results: Tapping the target with a pen


I have no ability in my lab to generate high speed projectiles..... The best I do now is to tap a target surface with a pen. But even this data is instructive....

Analysis Window: 164 μ s


- Equation 1 predicts a velocity resolution of .9 mm/s or 3 mm/s for 95% containment of the signal. This checks out!

I can get a 30 dB signal with about -30 dB light return.

- My sensitivity is about -50 dB (with a 10 dB signal requirement).




532 nm PDV data, 50 mW on target, -30 dB return




Next Step: Laser Driven Metal Foil

- MSTS LAO has agreed to let me conduct tests with their laser drive system
- Expect velocities of 500 m/s or greater with a 200 ps, 500 mJ pulse
- Can perform 10+ shots a day
- Will confirm Resolution and Sensitivity Performance (can compare with 1550 PDV).



Some Challenges

- Coupling of laser into 2.5 micron core not trivial.
- Greater attenuation in fiber (2-3 dB for a 50 m downlead, Nufern 460-HP).
- Less isolation (22 dB) and directivity (40 dB) and greater insertion loss (2.8 dB one-way) in Circulator (Ascentta, Inc.)
- Power density in core is high (~10 dB more than 1550 nm)
- **Sensitivity of detector (8 dB less A/W)**
- **BW of Detector (Commercially available detectors >10GHz are rare), limits max velocity to <3 km/s.**
- **Laser Classification: 3B at >5 mW.**



Conclusions

PDV at higher frequencies has several interesting attributes.

- Resolution of Motion, Sub-micron particle detection, Diffraction limited spot size, Higher critical density in plasmas.
- Niche applications are intriguing....

But there are challenges

- Perhaps sensitivity is the biggest concern – ~15 dB less sensitivity than 1550 nm PDV (could increase sensitivity by removing circulator).

Early benchtop results have demonstrated a useful sensitivity

- Next steps are experiments with a laser driven foil.
- Then repackaging and exploring applications for characterizing detonators, measuring particle size, and possibly studying WDM at DARHT.