

HiFiPDV:

A High Fidelity Analysis Program for Processing Single Velocity PDV Signals


Travis J. Voorhees
Sandia National Laboratories
Livermore, CA

PDV Workshop
February 7, 2023



Topics Covered (some will be quick)


- **Why was this code written?**
 - Who is this aimed towards?
 - Where will the code work well and not well?
- **How does the code work?**
- **Positives and negatives of this approach**
- **Areas of improvement for future releases**
- **Acknowledgements**
- **Q&A**



Why was this code written? Who is the audience?


Benefits of this code:

- **Semi-automated.** Minimal user inputs.
- **Consistent** analysis signal-to-signal.
- Calculates **systematic uncertainty.**
- **Simple to use.** Guides user at every step.




Who/What is this code for?

- People that are new to PDV analysis, i.e., graduate students, interns.
- Single velocity signals.
- >300 ns signals.
- Great for high noise:signal ratio.




There are a lot of knobs in PDV analysis codes

1. **Duration length**
 - a. Time-bin width, or number of data points in calculation
2. **Duration overlap**
 - a. Number of shared data points between adjacent durations
3. **Window function**
 - a. Tapering transform applied to DFT sine waves
4. **Peak extraction method**
 - a. Peak of DFT curve is the output velocity for a time-bin
5. **Smoothing**



Better **time** resolution




Better **velocity** resolution


There are a lot of knobs in PDV analysis codes

1. **Duration length**
 - a. Time-bin width, or number of data points in calculation
2. **Duration overlap**
 - a. Number of shared data points between adjacent durations
3. **Window function**
 - a. Tapering transform applied to DFT sine waves
4. **Peak extraction method**
 - a. Peak of DFT curve is the output velocity for a time-bin
5. **Smoothing**


Hamming window ($a_0 = 0.53836$)




Fourier transform



Hann window




Fourier transform




There are a lot of knobs in PDV analysis codes

1. **Duration length**
 - a. Time-bin width, or number of data points in calculation
2. **Duration overlap**
 - a. Number of shared data points between adjacent durations
3. **Window function**
 - a. Tapering transform applied to DFT sine waves
4. **Peak extraction method**
 - a. Peak of DFT curve is the output velocity for a time-bin
5. **Smoothing**





There are a lot of knobs in PDV analysis codes

1. **Duration length**
 - a. Time-bin width, or number of data points in calculation
2. **Duration overlap**
 - a. Number of shared data points between adjacent durations
3. **Window function**
 - a. Tapering transform applied to DFT sine waves
4. **Peak extraction method**
 - a. Peak of DFT curve is the output velocity for a time-bin
5. **Smoothing**



These knobs can significantly change output



**What knob settings are right for your signal?
What's the associated uncertainty?**


HiFiPDV's Solution: use a broad range of knob settings

Knob inputs:

1. Duration length
 - a. 5 to 100 ns
2. Duration overlap
 - a. 2.5 ns resolution
3. Window function
 - a. Hamming
 - b. Hann
 - c. Blackman
4. Peak Extraction method
 - a. Maximum
 - b. Robust centroid
 - c. Gaussian
5. Smoothing
 - a. Optional (output both)

Downsides to this approach:

1. Computational overhead




— How does *HiFiPDV* work? (mile high view)

Calculate Spectrograms

- 1 Time resolution (2.5 ns/pt)
- 1 Freq. resolution (2^{15} bins)
- 3 Window functions
- x20 Duration lengths
- 60 Spectrograms

User defines the ROI (white)

Example Spectrogram




Extract Velocities

- Max of spectral density
- Centroid of spectrum
- Gaussian fit to spectrum

180 velocities per time point


All Histories: Extraction Method



Evaluate Distribution Of Velocities

- Collapse histories into a heatmap
- Evaluate heatmap distribution at each time-bin for uncertainty analysis

Heatmap of All Histories



HiFiPDV tries to help guide the user. It has pop-up prompts at every step of the analysis

Data Import Settings

Time Units (s,us,ns):

Approximate Start Time (μ s):

Approximate End Time (μ s):


Laser Wavelength (nm):

Run Gaussian Fit? (yes/no)

Calculate 10-50-90 Rise Times (yes/no)?

Window Material:

Correct for frequency-shifting? (up/down/no)



```
Starting parallel pool (parpool) using the 'local' profile ...
Connected to the parallel pool (number of workers: 8).
User selected file: 22-4-024_S102C1.isf
Gaussian fitting process selected.
Signal is *NOT* frequency-shifted.
No window correction factor applied
Warning: All expected data was read from 22-4-024_S102C1.isf, but there still
> In isfread (line 165)
In HiFiPDV (line 302)
Number of data points per 1 ns is: 20.0
Point-to-point time resolution set to: 2.5000ns
Frequency bounds entered, beginning extraction process.
If the progress bar does not delete itself, please run the following commands:
First: >> F = findall(0,'type','figure','tag','TMWaitbar')
Then: >> delete(F)

=====
===== Extraction Process Finished! =====
=====

Total elapsed time: 2139.01 seconds
- Spectrogram production time: 27.96 seconds
- Maximum extraction method time: 2.57 seconds
- Robust Centroid extraction method time: 5.06 seconds
- Gaussian Peak extraction method time: 2103.21 seconds
```

Extracting Data from Signal

Extracting histories: Gaussian Method





Save Data?

Would you like to save your history, spectrogram, heatmap, and plots?

How well can *HiFiPDV* reproduce synthetic data?


Signal Parameters:

- Noise:Signal = 10
- Sample Rate = $4e10$
- Various accelerations
- Jagged or curved




Spectrograms of Signals


Synthetic Signal - A




Synthetic Signal - B



Synthetic Signal - C




Synthetic Signal - D




Final Code Output History


Synthetic Signal - A




Synthetic Signal - B







Synthetic Signal - C








Synthetic Signal - D



Overlaid Synthetic Input vs. Code Output



Applied to real data: Cu @ 2.39 km/s into Mo-PFO RMI



Areas of improvement


- Computational overhead
- Speed up Gaussian fit
 - Gaussian fitting accounts for ~95-99% of code run time
- Better spectrogram cropping
- More output freedom
 - Choose your spectrogram
 - Choose your colormaps
- Reduce toolbox usage



```
=====  
===== Extraction Process Finished! =====  
=====  
Total elapsed time: 2139.01 seconds  
- Spectrogram production time: 27.96 seconds  
- Maximum extraction method time: 2.57 seconds  
- Robust Centroid extraction method time: 5.06 seconds  
- Gaussian Peak extraction method time: 2103.21 seconds
```

Acknowledgements

- Code Testers:
 - Georgia Tech: Greg Kennedy, Andrew Boddorff, Katie Koube, Keara Frawley, Ben Zussman, and the whole HSR lab
 - Sandia Natl Lab: Stacy Guo
- Advice and Support:
 - Sandia Natl Lab: Dan Dolan, Sheri Payne
 - Georgia Tech: Greg Kennedy, Naresh Thadhani
 - Los Alamos Natl Lab: John Lang
- Example Data:
 - Georgia Tech HSR Lab
 - Dynamic Compression Sector (DCS)



Thank you!
tjvoorh@sandia.gov

Inspired by Dan Dolan's "Accuracy and Precision in PDV" Influenced by SIRHEN, PDV_Display, and PlotData


Dan's methodology: 10.1063/1.3429257

The following procedure is used to simulate the resolution of single-signal PDV measurements.


- (1) A random phase δ is selected from the range $(0, 2\pi)$.
- (2) The sampled signal s_k is calculated from Eq. (3) for the current value of f_0 and σ (drawing R_k from a normal distribution centered at zero with a standard deviation of unity).
- (3) The power spectrum $P(f)$ is determined from a fast Fourier transform (FFT) (Ref. 7) of s_k using boxcar/rectangle, Hamming, and Hann window functions.⁸
- (4) The peak location \bar{f} is determined by fitting $P(f)$ with a Gaussian curve.
- (5) Steps 1–4 are repeated 20 000 times to create a distribution of \bar{f} values for each window function at the current value of f and σ .
- (6) Accuracy is calculated from the difference between $\langle \bar{f} \rangle$ and f_0 .
- (7) Precision is calculated from the standard deviation of \bar{f} .

How does the code work?

1. Import and rough crop
2. Spectrograms & Cropping:
 - a. Generates spectrograms
 - b. Users define ROI
 - c. Crops specs to ROI
3. Extracts histories from cropped spectrograms
4. Collapses all histories into a single heatmap
5. Extracts final history and 95% confidence interval from the heatmap




How are duration length and overlap ratio coupled?








- At each time location, vary:
- Duration length, DN
 - Tapering window function
 - Hamming
 - Hann
 - Blackman






How does the code work?








Shot1907, Probe 1: Copper-Baseplate



Shot1930, Probe MockAMEM: Example




Shot1707, Probe TiB_2 -Confinement: Low-SNR-Example




Shot APS-19-1-007, Probe S100C1: Multiple-Velocity-Example


Example Spectrogram




Final Output History and Uncertainties



All Histories: Extraction Method



All Histories: Window Function



Heatmap of All Histories

